

		List of Compatibility for ANKYLOS ^(R) C/X	
		ENGLISH	ENGLISH
SAP	BaaN	Product description	compatible with ANKYLOS C/X
31010080	17-1010	Cover Screw (spare part)	not compatible with ANYKLOS C/X
31010110	17-1040	C/X Membrane Screw Ø 3,5 mm	compatible with ANKYLOS C/X
31010120	17-1060	C/X Membrane Screw Ø 6,0 mm	compatible with ANKYLOS C/X
31010130	17-1061	C/X Membrane Screw Ø 6,0 mm for sinus lift, 1mm distance	compatible with ANKYLOS C/X
31010140	17-1062	C/X Membrane Screw Ø 6,0 mm for sinus lift, 2mm distance	compatible with ANKYLOS C/X, new product
31010405	17-0542	ANKYLOS C/X Impl A8/ø3,5/L8	new product for ANKYLOS C/X
31010408	17-0543	ANKYLOS C/X ImplA9,5/ø3,5/L9,5	new product for ANKYLOS C/X
31010410	17-0544	ANKYLOS C/X Impl A11/ø3,5/L11	new product for ANKYLOS C/X
31010420	17-0546	ANKYLOS C/X Impl A14/ø3,5/L14	new product for ANKYLOS C/X
31010422	17-0547	ANKYLOS C/X Impl A17/ø3,5/L17	new product for ANKYLOS C/X
31010425	17-0552	ANKYLOS C/X Impl B8/ø4,5/L8	new product for ANKYLOS C/X
31010428	17-0553	ANKYLOS C/X ImplB9,5/ø4,5/L9,5	new product for ANKYLOS C/X
31010430	17-0554	ANKYLOS C/X Impl B11/ø4,5/L11	new product for ANKYLOS C/X
31010440	17-0556	ANKYLOS C/X Impl B14/ø4,5/L14	new product for ANKYLOS C/X
31010450	17-0557	ANKYLOS C/X Impl B17/ø4,5/L17	new product for ANKYLOS C/X
31010453	17-0562	ANKYLOS C/X Impl C8/ø5,5/L8	new product for ANKYLOS C/X
31010455	17-0563	ANKYLOS C/X ImplC9,5/ø5,5/L9,5	new product for ANKYLOS C/X
31010458	17-0564	ANKYLOS C/X Impl C11/ø5,5/L11	new product for ANKYLOS C/X
31010460	17-0566	ANKYLOS C/X Impl C14/ø5,5/L14	new product for ANKYLOS C/X
31010470	17-0567	ANKYLOS C/X Impl C17/ø5,5/L17	new product for ANKYLOS C/X
31010472	17-0572	ANKYLOS C/X Impl D8/ø7,0/L8	new product for ANKYLOS C/X
31010474	17-0573	ANKYLOS C/X ImplD9,5/ø7,0/L9,5	new product for ANKYLOS C/X
31010476	17-0574	ANKYLOS C/X Impl D11/ø7,0/L11	new product for ANKYLOS C/X
31010478	17-0576	ANKYLOS C/X Impl D14/ø7,0/L14	new product for ANKYLOS C/X
31010480	17-1001	ANKYLOS C/X Cover screw	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31021010	17-1131	Sulcus Former for Standard abutments a/1,5	not compatible with ANYKLOS C/X
31021011	17-1132	Standard C/ Sulcus Former compact for abutment a/1,5	only compatible, when laser marked with "C/"
31021020	17-1133	Sulcus Former for Standard abutments a/3,0	not compatible with ANYKLOS C/X
31021021	17-1134	Standard C/ Sulcus Former compact for abutment a/3,0	only compatible, when laser marked with "C/"
31021026	17-1136	Standard C/ Sulcus Former compact for abutment a/4,5	only compatible, when laser marked with "C/"
31021028	17-1138	Standard C/ Sulcus Former compact for abutment a/6,0	only compatible, when laser marked with "C/"
31021030	17-1151	Sulcus Former for Standard abutments b/1,5	not compatible with ANYKLOS C/X
31021031	17-1152	Standard C/ Sulcus Former compact for abutment b/1,5	only compatible, when laser marked with "C/"

List of Compatibility for ANKYLOS ^(R) C/X			
		ENGLISH	ENGLISH
SAP	BaaN	Product description	compatible with ANKYLOS C/X
31021060	17-2132	Standard C/ Abutment a/1,5/6,0 straight	only compatible, when laser marked with "C/"
31021070	17-2133	Standard C/ Abutment a/3,0/4,0 straight	only compatible, when laser marked with "C/"
31021076	17-2134	Standard C/ Abutment a/4,5/4,0 straight	only compatible, when laser marked with "C/"
31021078	17-2135	Standard C/ Abutment a/6,0/4,0 straight	only compatible, when laser marked with "C/"
31021080	17-2136	Standard C/ Abutment a/3,0/6,0 straight	only compatible, when laser marked with "C/"
31021090	17-2151	Standard C/ Abutment b/1,5/4,0 straight	only compatible, when laser marked with "C/"
31021100	17-2152	Standard C/ Abutment b/1,5/6,0 straight	only compatible, when laser marked with "C/"
31021110	17-2153	Standard C/ Abutment b/3,0/4,0 straight	only compatible, when laser marked with "C/"
31021120	17-2156	Standard C/ Abutment b/3,0/6,0 straight	only compatible, when laser marked with "C/"
31021130	17-2231	Standard C/ Abutment a/1,5/4,0 15°	only compatible, when laser marked with "C/"
31021140	17-2232	Standard C/ Abutment a/1,5/6,0 15°	only compatible, when laser marked with "C/"
31021150	17-2233	Standard C/ Abutment a/3,0/4,0 15°	only compatible, when laser marked with "C/"
31021160	17-2236	Standard C/ Abutment a/3,0/6,0 15°	only compatible, when laser marked with "C/"
31021170	17-2251	Standard C/ Abutment b/1,5/4,0 15°	only compatible, when laser marked with "C/"
31021180	17-2252	Standard C/ Abutment b/1,5/6,0 15°	only compatible, when laser marked with "C/"
31021190	17-2253	Standard C/ Abutment b/3,0/4,0 15°	only compatible, when laser marked with "C/"
31021200	17-2256	Standard C/ Abutment b/3,0/6,0 15°	only compatible, when laser marked with "C/"
31021300	17-2331	Snap Attachment C/ a/1,5 complete	only compatible, when laser marked with "C/"
31021301	17-2330	Snap Attachment C/ a/1,5	only compatible, when laser marked with "C/"
31021310	17-2333	Snap Attachment C/ a/3,0 complete	only compatible, when laser marked with "C/"
31021311	17-2332	Snap Attachment C/ a/3,0	only compatible, when laser marked with "C/"
31021312	17-2339	Snap for snap attachment	compatible, no connection to implant
31021313	17-2338	Silicone ring for snap attachment	compatible, no connection to implant
31021320	17-2335	Snap Attachment C/ a/4,5 complete	only compatible, when laser marked with "C/"
31021321	17-2334	Snap Attachment C/ a/4,5	only compatible, when laser marked with "C/"
31021400	17-2451	Abutment with Magnetic Attachment b/1,5 complete	not compatible with ANYKLOS C/X, product terminated
31021401	17-2450	Abutment with Magnetic Attachment b/1,5	not compatible with ANYKLOS C/X, product terminated
31021402	17-2459	Secondary magnet	product terminated
31021403	17-2458	Positioning sleeve supra for secondary magnet	not compatible with ANYKLOS C/X, product terminated
31021404	17-2457	Positioning sleeve sub for secondary magnet	not compatible with ANYKLOS C/X, product terminated
31021405	17-2750	Polymerisation sleeve for SynCone	compatible, no connection to implant
31021410	17-2453	Abutment with Magnetic Attachment b/3,0 complete	product terminated
31021411	17-2452	Abutment with Magnetic Attachment b/3,0	product terminated

		List of Compatibility for ANKYLOS ^(R) C/X	
		ENGLISH	ENGLISH
SAP	BaaN	Product description	compatible with ANKYLOS C/X
31021450	17-1941	Sulcus former A for Permador abutment	not compatible with ANKYLOS C/X
31021451	17-1942	Sulcus former A for Permador abutment C/	new product for ANKYLOS C/X, not compatible with ANKYLOS plus
31021452	17-1951	Sulcus former B for Permador abutment	compatible with ANKYLOS C/X
31021454	17-1961	Sulcus former C for Permador abutment	compatible with ANKYLOS C/X
31021505	17-4111	Screw for Sulcus former Balance Anterior 1,5	compatible with ANKYLOS C/X
31021506	17-4113	Screw for temporary crown Balance Anterior 1,5	product terminated
31021507	17-4112	Screw for Sulcus former Balance Anterior 3,0	compatible with ANKYLOS C/X
31021508	17-4114	Screw for temporary crown Balance Anterior 3,0	product terminated
31021510	17-1262	Balance Anterior C/ Sulcus former 1,5 small (with screw)	compatible with ANKYLOS C/X
31021515	17-1263	Balance Anterior C/ Sulcus former 3,0 small (with screw)	compatible with ANKYLOS C/X
31021520	17-1272	Balance Anterior C/ Sulcus former 1,5 large (with screw)	compatible with ANKYLOS C/X
31021525	17-1273	Balance Anterior C/ Sulcus former 3,0 large (with screw)	compatible with ANKYLOS C/X
31021528	17-1251	Balance Posterior C/ Sulcus former 0,75	only compatible, when laser marked with "C/"
31021530	17-1252	Balance Posterior C/ Sulcus former 1,5	only compatible, when laser marked with "C/"
31021535	17-1253	Balance Posterior C/ Sulcus former 3,0	only compatible, when laser marked with "C/"
31021540	17-1254	Balance Posterior C/ Sulcus former 4,5	only compatible, when laser marked with "C/"
31021550	17-3161	Balance Anterior abutment C/ small / 1,5 straight	compatible with ANKYLOS C/X
31021560	17-3162	Balance Anterior abutment C/ small / 3,0 straight	compatible with ANKYLOS C/X
31021570	17-3171	Balance Anterior abutment C/ large / 1,5 straight	compatible with ANKYLOS C/X
31021580	17-3172	Balance Anterior abutment C/ large / 3,0 straight	compatible with ANKYLOS C/X
31021590	17-3261	Balance Anterior abutment C/ small / 1,5 15°	compatible with ANKYLOS C/X
31021600	17-3262	Balance Anterior abutment C/ small / 3,0 15°	compatible with ANKYLOS C/X
31021610	17-3271	Balance Anterior abutment C/ large / 1,5 15°	compatible with ANKYLOS C/X
31021620	17-3272	Balance Anterior abutment C/ large / 3,0 15°	compatible with ANKYLOS C/X
31021625	17-3151	Balance Posterior abutment C/ 0,75 straight	compatible with ANKYLOS C/X
31021630	17-3152	Balance Posterior abutment C/ 1,5 straight	compatible with ANKYLOS C/X
31021640	17-3153	Balance Posterior abutment C/ 3,0 straight	compatible with ANKYLOS C/X
31021645	17-3154	Balance Posterior abutment C/ 4,5 straight	compatible with ANKYLOS C/X
31021647	17-3251	Balance Posterior abutment C/ 0,75 7,5°	compatible with ANKYLOS C/X
31021648	17-3351	Balance Posterior abutment C/ 0,75 15°	compatible with ANKYLOS C/X
31021649	17-3252	Balance Posterior abutment C/ 1,5 7,5°	compatible with ANKYLOS C/X
31021650	17-3352	Balance Posterior abutment C/ 1,5 15°	compatible with ANKYLOS C/X
31021652	17-3451	Balance Posterior abutment C/ 0,75 22,5°	compatible with ANKYLOS C/X

List of Compatibility for ANKYLOS ^(R) C/X			
		ENGLISH	ENGLISH
SAP	BaaN	Product description	compatible with ANKYLOS C/X
31021659	17-3651	Balance Posterior abutment C/ 0,75 37,5°	compatible with ANKYLOS C/X
31021660	17-3353	Balance Posterior abutment C/ 3,0 15°	compatible with ANKYLOS C/X
31021665	17-3354	Balance Posterior abutment C/ 4,5 15°	compatible with ANKYLOS C/X
31021668	17-3751	Balance Base Abutment C/ 0,75	only compatible, when laser marked with "C/"
31021670	17-3752	Balance Base Abutment C/ 1,5	only compatible, when laser marked with "C/"
31021680	17-3753	Balance Base Abutment C/ 3,0	only compatible, when laser marked with "C/"
31021690	17-3754	Balance Base Abutment C/ 4,5	only compatible, when laser marked with "C/"
31021710	17-2741	SynCone C/ Abutment 6° conical head 1,5 / 0°	compatible with ANKYLOS C/X
31021712	17-2541	SynCone C/ Abutment 4° conical head 1,5 / 0°	compatible with ANKYLOS C/X
31021720	17-2742	SynCone C/ Abutment 6° conical head 3,0 / 0°	compatible with ANKYLOS C/X
31021722	17-2542	SynCone C/ Abutment 4° conical head 3,0 / 0°	compatible with ANKYLOS C/X
31021730	17-2743	SynCone C/ Abutment 6° conical head 4,5 / 0°	compatible with ANKYLOS C/X
31021732	17-2543	SynCone C/ Abutment 4° conical head 4,5 / 0°	compatible with ANKYLOS C/X
31021750	17-2748	SynCone Cap 6° conical head with retention	compatible, no connection to implant
31021752	17-2548	SynCone Cap 4° conical head with retention	compatible, no connection to implant
31021770	17-2841	SynCone C/ Abutment 6° conical head 1,5/ 15°	compatible with ANKYLOS C/X
31021772	17-2641	SynCone C/ Abutment 4° conical head 1,5/ 15°	compatible with ANKYLOS C/X
31021780	17-2842	SynCone C/ Abutment 6° conical head 3,0/ 15°	compatible with ANKYLOS C/X
31021782	17-2642	SynCone C/ Abutment 4° conical head 3,0/ 15°	compatible with ANKYLOS C/X
31021790	17-2843	SynCone C/ Abutment 6° conical head 4,5/ 15°	compatible with ANKYLOS C/X
31021792	17-2643	SynCone C/ Abutment 4° conical head 4,5/ 15°	compatible with ANKYLOS C/X
31021810	17-3361	Cercon Balance small 1,5 straight	compatible with ANKYLOS C/X
31021820	17-3362	Cercon Balance small 3,0 straight	compatible with ANKYLOS C/X
31021830	17-3371	Cercon Balance large 1,5 straight	compatible with ANKYLOS C/X
31021840	17-3372	Cercon Balance large 3,0 straight	compatible with ANKYLOS C/X
31021850	17-3461	Cercon Balance small 1,5 15°	compatible with ANKYLOS C/X
31021860	17-3462	Cercon Balance small 3,0 15°	compatible with ANKYLOS C/X
31021870	17-3471	Cercon Balance large 1,5 15°	compatible with ANKYLOS C/X
31021880	17-3472	Cercon Balance large 3,0 15°	compatible with ANKYLOS C/X
31021912	17-2941	SynCone C/ Abutment 4° conical head 1,5/ 22,5°	compatible with ANKYLOS C/X
31023375	17-7073	Cercon Balance Anterior abutment C/ large 3.0 15°, dentin	compatible with ANKYLOS C/X
31023370	17-7072	Cercon Balance Anterior abutment C/ large 1.5 15°, dentin	compatible with ANKYLOS C/X
31023355	17-7069	Cercon Balance Anterior abutment C/ large 3.0 straight, dentin	compatible with ANKYLOS C/X

List of Compatibility for ANKYLOS ^(R) C/X			
		ENGLISH	ENGLISH
SAP	BaaN	Product description	compatible with ANKYLOS C/X
31023315	17-7061	Cercon Balance Anterior abutment C/ small 3.0 straight, dentin	compatible with ANKYLOS C/X
31023310	17-7060	Cercon Balance Anterior abutment C/ small 1.5 straight, dentin	compatible with ANKYLOS C/X
31023075	17-7014	Cercon Balance Anterior abutment C/ large 3.0 15°, neutral	compatible with ANKYLOS C/X
31023070	17-7013	Cercon Balance Anterior abutment C/ large 1.5 15°, neutral	compatible with ANKYLOS C/X
31023055	17-7010	Cercon Balance Anterior abutment C/ large 3.0 straight, neutral	compatible with ANKYLOS C/X
31023050	17-7009	Cercon Balance Anterior abutment C/ large 1.5 straight, neutral	compatible with ANKYLOS C/X
31023035	17-7006	Cercon Balance Anterior abutment C/ small 3.0 15°, neutral	compatible with ANKYLOS C/X
31023030	17-7005	Cercon Balance Anterior abutment C/ small 1.5 15°, neutral	compatible with ANKYLOS C/X
31023015	17-7002	Cercon Balance Anterior abutment C/ small 3.0 straight, neutral	compatible with ANKYLOS C/X
31023010	17-7001	Cercon Balance Anterior abutment C/ small 1.5 straight, neutral	compatible with ANKYLOS C/X
31023336	17-7265	Balance C/ Temporary abutment small 3.0 15° dentin	compatible with ANKYLOS C/X
31023316	17-7261	Balance C/ Temporary abutment small 3.0 0° dentin	compatible with ANKYLOS C/X
31023331	17-7264	Balance C/ Temporary abutment small 1.5 15° dentin	compatible with ANKYLOS C/X
31023311	17-7260	Balance C/ Temporary abutment small 1.5 0° dentin	compatible with ANKYLOS C/X
31023376	17-7273	Balance C/ Temporary abutment large 3.0 15° dentin	compatible with ANKYLOS C/X
31023356	17-7269	Balance C/ Temporary abutment large 3.0 0° dentin	compatible with ANKYLOS C/X
31023371	17-7272	Balance C/ Temporary abutment large 1.5 15° dentin	compatible with ANKYLOS C/X
31023351	17-7268	Balance C/ Temporary abutment large 1.5 0° dentin	compatible with ANKYLOS C/X
31024010	17-1301	Regular C/X Gingiva former GH 0,75	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024012	17-1302	Regular C/X Gingiva former GH 1,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024014	17-1303	Regular C/X Gingiva former GH 3,0	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024016	17-1304	Regular C/X Gingiva former GH 4,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024110	17-3010	Regular /X Abutment GH 0,75 A 0	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024111	17-3011	Regular /X Abutment GH 0,75 A 7,5	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024112	17-3012	Regular /X Abutment GH 0,75 A 15	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024113	17-3013	Regular /X Abutment GH 0,75 A 22,5	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024114	17-3014	Regular /X Abutment GH 0,75 A 30	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024115	17-3015	Regular /X Abutment GH 0,75 A 37,5	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024120	17-3020	Regular /X Abutment GH 1,5 A 0	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024121	17-3021	Regular /X Abutment GH 1,5 A 7,5	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024122	17-3022	Regular /X Abutment GH 1,5 A 15	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024123	17-3023	Regular /X Abutment GH 1,5 A 22,5	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024124	17-3024	Regular /X Abutment GH 1,5 A 30	new product for ANKYLOS C/X, not compatible with ANYKLOS plus

		List of Compatibility for ANKYLOS ^(R) C/X	
		ENGLISH	ENGLISH
SAP	BaaN	Product description	compatible with ANKYLOS C/X
31024132	17-3032	Regular /X Abutment GH 3,0 A 15	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024133	17-3033	Regular /X Abutment GH 3,0 A 22,5	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024134	17-3034	Regular /X Abutment GH 3,0 A 30	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024135	17-3035	Regular /X Abutment GH 3,0 A 37,5	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024140	17-3040	Regular /X Abutment GH 4,5 A 0	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024141	17-3041	Regular /X Abutment GH 4,5 A 7,5	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024142	17-3042	Regular /X Abutment GH 4,5 A 15	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024143	17-3043	Regular /X Abutment GH 4,5 A 22,5	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024144	17-3044	Regular /X Abutment GH 4,5 A 30	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024145	17-3045	Regular /X Abutment GH 4,5 A 37,5	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31024150	17-3050	Regular C/ Abutment GH 0,75 A 0	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024151	17-3051	Regular C/ Abutment GH 0,75 A 7,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024152	17-3052	Regular C/ Abutment GH 0,75 A 15	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024153	17-3053	Regular C/ Abutment GH 0,75 A 22,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024154	17-3054	Regular C/ Abutment GH 0,75 A 30	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024155	17-3055	Regular C/ Abutment GH 0,75 A 37,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024160	17-3060	Regular C/ Abutment GH 1,5 A 0	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024161	17-3061	Regular C/ Abutment GH 1,5 A 7,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024162	17-3062	Regular C/ Abutment GH 1,5 A 15	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024163	17-3063	Regular C/ Abutment GH 1,5 A 22,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024164	17-3064	Regular C/ Abutment GH 1,5 A 30	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024165	17-3065	Regular C/ Abutment GH 1,5 A 37,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024170	17-3070	Regular C/ Abutment GH 3,0 A 0	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024171	17-3071	Regular C/ Abutment GH 3,0 A 7,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024172	17-3072	Regular C/ Abutment GH 3,0 A 15	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024173	17-3073	Regular C/ Abutment GH 3,0 A 22,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024174	17-3074	Regular C/ Abutment GH 3,0 A 30	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024175	17-3075	Regular C/ Abutment GH 3,0 A 37,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024180	17-3080	Regular C/ Abutment GH 4,5 A 0	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024181	17-3081	Regular C/ Abutment GH 4,5 A 7,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024182	17-3082	Regular C/ Abutment GH 4,5 A 15	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024183	17-3083	Regular C/ Abutment GH 4,5 A 22,5	new product for ANKYLOS C/X, compatible with ANYKLOS plus
31024184	17-3084	Regular C/ Abutment GH 4,5 A 30	new product for ANKYLOS C/X, compatible with ANYKLOS plus

List of Compatibility for ANKYLOS ^(R) C/X			
		ENGLISH	ENGLISH
SAP	BaaN	Product description	compatible with ANKYLOS C/X
31033093	17-5056	Insertion Instrument for implants, long	not compatible, new insertion post
31033101	17-5051	Insertion Instrument for implants, short	not compatible, new insertion post
31033102	17-5061	Insertion Instrument for implants, short (M)	not compatible, new insertion post
31033110	17-5042	Ratchet-insert for implants, medium	not compatible, new insertion post
31033111	17-5052	Insertion Instrument for implants, medium	product terminated
31033120	17-5043	Ratchet insert for implants, long	not compatible, new insertion post
31033121	17-5053	Insertion Instrument for implants, long	not compatible, new insertion post
31036000	17-8010	Surgery Tray for A-Implants	compatible, if equipped with new instruments
31036010	17-8020	Surgery Tray for A+B-Implants	compatible, if equipped with new instruments
31036011	17-8021	Surgery Tray for A+B-Implants (M)	compatible, if equipped with new instruments
31036020	17-8030	Surgery Tray for A+B+C-Implants	compatible, if equipped with new instruments
31036025	17-8040	Surgery Tray for D-Implants	compatible, if equipped with new instruments
31036500	17-8200	Surgical Kit manual ABC	compatible, if equipped with new instruments
31036501	17-8201	ANKYLOS C/X Surgical Kit ABC	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31036510	17-8210	Surgical Kit manual AB	compatible, if equipped with new instruments
31036511	17-8211	ANKYLOS C/X Surgical Kit AB	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31036521	17-8221	ANKYLOS C/X Surgical Kit motor AB	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31036522	17-8222	Surgical Kit motor AB reamer	compatible, if equipped with new instruments
31037502	17-5302	ANKYLOS C/X Open-end wrench	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31037510	17-5014	Ratchet insert for M-instruments short HXSS	compatible, no relevance for index
31037511	17-5012	Ratchet insert for M-instruments medium HXSS	compatible, no relevance for index
31037512	17-5013	Ratchet insert for M-instruments long HXSS	compatible, no relevance for index
31037531	17-5064	Insertion instrument for C/X implants, manual short	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31037532	17-5065	Insertion instrument for C/X implants, manual medium	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31037533	17-5066	Insertion instrument for C/X implants, manual long	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31037539	17-5058	Insertion handle for C/X implants, extra long	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31037541	17-5067	Insertion instrument for C/X implants HXSS motor, short	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31037543	17-5069	Insertion instrument for C/X implants HXSS motor, long	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31045270	17-4301	Balance C/ Implant analog A	only compatible, when restoration without index
31045272	17-4302	Balance C/ Implant analog B	only compatible, when restoration without index
31045274	17-4303	Balance C/ Implant analog C	only compatible, when restoration without index
31045400	17-1701	Transfer screw, short	only compatible, when restoration without index
31045402	17-1702	Transfer screw, long	only compatible, when restoration without index

		List of Compatibility for ANKYLOS ^(R) C/X	
		ENGLISH	ENGLISH
SAP	BaaN	Product description	compatible with ANKYLOS C/X
31045416	17-1714	Repositioning post for Balance abutments, narrow	only compatible, when laser marked with "C/"
31045710	17-8565	Laboratory Trial Abutment Tray Balance Posterior	compatible, no relevance for index
31047011	17-4221	Laboratory implant C/X A	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31047012	17-4222	Laboratory implant C/X B	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31047013	17-4223	Laboratory implant C/X C	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31047020	17-1800	C/X Transfer post PickUp technique, short	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31047023	17-1802	C/X Transfer post PickUp technique, long	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31047022	17-1801	PickUp screw extension	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31047030	17-1806	C/X Reposition post	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31047110	17-1810	Regular C/X Abutment analog A0	new product for ANKYLOS C/X, to be used with Regular abutments only
31047111	17-1811	Regular C/X Abutment analog A7,5	new product for ANKYLOS C/X, to be used with Regular abutments only
31047112	17-1812	Regular C/X Abutment analog A15	new product for ANKYLOS C/X, to be used with Regular abutments only
31047113	17-1813	Regular C/X Abutment analog A22,5	new product for ANKYLOS C/X, to be used with Regular abutments only
31047114	17-1814	Regular C/X Abutment analog A30	new product for ANKYLOS C/X, to be used with Regular abutments only
31047115	17-1815	Regular C/X Abutment analog A37,5	new product for ANKYLOS C/X, to be used with Regular abutments only
31047120	17-4230	Regular C/X 3-in-1 coping A0	new product for ANKYLOS C/X, to be used with Regular abutments only
31047121	17-4231	Regular C/X 3-in-1 coping A7,5	new product for ANKYLOS C/X, to be used with Regular abutments only
31047122	17-4232	Regular C/X 3-in-1 coping A15	new product for ANKYLOS C/X, to be used with Regular abutments only
31047123	17-4233	Regular C/X 3-in-1 coping A22,5	new product for ANKYLOS C/X, to be used with Regular abutments only
31047124	17-4234	Regular C/X 3-in-1 coping A30	new product for ANKYLOS C/X, to be used with Regular abutments only
31047125	17-4235	Regular C/X 3-in-1 coping A37,5	new product for ANKYLOS C/X, to be used with Regular abutments only
31056010	17-4131	Headless Screw	not compatible with Regular abutment
31056021	17-4132	Retaining Screw, occlusal hexagon M1,6	not compatible with Regular abutment
31056051	17-4133	Retaining Screw, lateral hexagon M1x0,2	not compatible with Regular Abutment
31056052	17-4134	Sleeve for retaining screw, lateral hexagon	not compatible with Regular abutment
31056054	17-4135	Retaining Screw lateral hexagon M1,2	not compatible with Regular abutment
31056100	17-2511	Gold Coping, Permador a/1,5/4,0 & a/3,0/4,0 with screw	compatible, no connection to implant
31056110	17-2512	Gold Coping, Permador a/1,5/6,0 & a/3,0/6,0 with screw	compatible, no connection to implant
31056120	17-2513	Gold Coping, Permador b/1,5/4,0 & b/3,0/4,0 with screw	compatible, no connection to implant
31056130	17-2514	Gold Coping, Permador b/1,5/6,0 & b/3,0/6,0 with screw	compatible, no connection to implant
31056140	17-4136	Screw for Gold Coping	compatible, no connection to implant
31056160	17-2521	Gold Coping, Degunorm a/1,5/4,0 & a/3,0/4,0 with screw	compatible, no connection to implant
31056170	17-2522	Gold Coping, Degunorm a/1,5/6,0 & a/3,0/6,0 with screw	compatible, no connection to implant

		List of Compatibility for ANKYLOS ^(R) C/X	
		ENGLISH	ENGLISH
SAP	BaaN	Product description	compatible with ANKYLOS C/X
31056210	17-3541	Gold Coping f. Balance Base Abut. (Degunorm) with screw	compatible, no connection to implant
31056215	17-3543	Retention Coping f. Balance Base Abutment	compatible, no connection to implant
31056220	17-2611	Threaded Sleeve for individual screw retention, short	compatible, no connection to implant
31056230	17-2612	Threaded sleeve for individual screw retention, medium	compatible, no connection to implant
31056240	17-2613	Threaded sleeve for individual screw retention, long	compatible, no connection to implant
31056250	17-2621	Splinting Screw for individual screw retention, short	compatible, no connection to implant
31056260	17-2622	Splinting screw for individual screw retention, medium	compatible, no connection to implant
31056270	17-2623	Splinting screw for individual screw retention, long	compatible, no connection to implant
31056301	17-4137	Fastening screw lateral hexagon M1,4	new product for ANKYLOS C/X, only for Regular abutments, not compatible with Balance abutments
31056304	17-4138	Horizontal screw seat for fixation screw M1,4	new product for ANKYLOS C/X, only for Regular abutments, not compatible with Balance abutments
31099011	17-9945	ANKYLOS C/X training implant A11	new product for ANKYLOS C/X, not compatible with ANYKLOS plus
31099200	17-9910	Acrylic model maxilla, crown & bridge, 3 impl.	not compatible with ANYKLOS C/X
31099240	17-9911	Acrylic model mandible edentulous, 4 impl.	not compatible with ANYKLOS C/X
31099260	17-9912	Acrylic model mandible edentulous, SynCone	not compatible with ANYKLOS C/X
31099261	17-9913	Acrylic model mandible edentulous, SynCone, sales mod.	not compatible with ANYKLOS C/X
31099262	17-9914	Acrylic model mandible SynCone, alignment, sales mod.	not compatible with ANYKLOS C/X
31099351	17-9961	Model 5:1 Ankylos C/X Implant A 11	new product for ANKYLOS C/X
31099352	17-9962	Model 5:1 Cover screw ANKYLOS C/X	new product for ANKYLOS C/X
31099353	17-9963	Model 5:1 ANKYLOS Regular C/X Gingiva former 3,0	new product for ANKYLOS C/X
31099354	17-9964	Model 5:1 ANKYLOS Regular _/X Abutment 3,0 A15	new product for ANKYLOS C/X
31099355	17-9965	Model 5:1 ANKYLOS Regular C/_ Abutment 3,0 A15	new product for ANKYLOS C/X
31099356	17-9966	Model 5:1 ANKYLOS Transfer abutment C/X, short	new product for ANKYLOS C/X
31099357	17-9967	Model 5:1 screw driver	new product for ANKYLOS C/X